

## SEQUENTIAL PROCESSING REACTION VESSEL FOR CHEMICAL FRACTIONATION AND ANALYSIS

A sequential processing reactor vessel and method is disclosed for accelerated  
5 extraction and fractionation of chemical analytes from complex solid sample materials. The  
device and method provide for sequential extraction of elemental constituents from solid  
materials by sequentially contacting target samples within a single reaction vessel using a  
series of different reagents at temperatures as high as 150° C and pressures up to 150 psi to  
accelerate reactions. The aggressive chemical treatments provided by the disclosed device  
10 and method provide for complete digestion of solid samples in liquid analyte samples that  
can be directly analyzed by conventional spectrometry or other suitable methods. The  
disclosed device and method provide for efficient sample processing and accelerated  
reactions to significantly reduce processing times and costs for elemental analysis of solids  
while improving accuracy, precision and reliability of results compared to conventional  
15 devices and methods. The disclosed device and method are compatible with both  
conventional convection and radiant heating sources as well as microwave heating and can  
be readily adapted to marine, geological, environmental, industrial and research solids  
analysis applications.